




Histone H3K14me1 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA923350 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | GeneID:8290SwissProt:Q16695 |
| Immunogen | A synthetic peptide of Human H3K14me1 |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB,IHC,IF;WB:1:500-1:2000,IHC:1:50-1:200,IF:1:20-1:50 |
| Form | Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3. |
| Purification Method | Antigen Affinity Purified |
| Clonality | Polyclonal |
| Alias | H3t; H3.4; H3/g; H3FT; |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | HIST3H3 |

Image



Western blot analysis of extracts of HeLa and H3 protein, using H3K14me1 antibody.

| | |
|-----------------------|---|
| Product Modify | K14me1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |