




Histone H3K14me2 Antibody

Product Code	CSB-PA830106
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	GeneID:8290SwissProt:Q16695
Immunogen	A synthetic peptide of Human H3K14me2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:2000,IHC:1:50-1:100,IF:1:20-1:50
Form	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	H3t; H3.4; H3/g; H3FT;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HIST3H3

Image



Western blot analysis of extracts of HeLa and H3 protein, using H3K14me2 antibody.

Product Modify	K14me2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.