




# Histone H3K14me3 Antibody

<b>Product Code</b>	CSB-PA256630
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	GeneID:8290SwissProt:Q16695
<b>Immunogen</b>	A synthetic peptide of Human H3K14me3
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB,IHC,IF;WB:1:500-1:2000,IHC:1:50-1:200,IF:1:20-1:50
<b>Form</b>	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Clonality</b>	Polyclonal
<b>Alias</b>	H3t; H3.4; H3/g; H3FT;
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	HIST3H3

## Image



Western blot analysis of extracts of HeLa and H3 protein, using H3K14me3 antibody.

<b>Product Modify</b>	K14me3
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.