




FOXO1 Antibody

Product Code	CSB-PA168688
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12778
Immunogen	A synthesized peptide derived from Human FOXO1A
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	FOXO1A Antibody detects endogenous levels of total FOXO1A
Tested Applications	ELISA,WB,IHC;WB:1:500-1:2000,IHC:1:50-1:200
Form	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	FOXO1A, FOXO1, Forkhead box protein O1A, Forkhead in rhabdomyosarcoma
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FOXO1

Image



Western blot analysis of FOXO1A expression in Serum treated HeLa whole cell lysates, The lane on the left is treated with the antigen-specific peptide