




# GSK3A Antibody

<b>Product Code</b>	CSB-PA205981
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P49840
<b>Immunogen</b>	A synthesized peptide derived from Human GSK3 alpha/beta
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Specificity</b>	GSK3 alpha/beta Antibody detects endogenous levels of total GSK3 alpha/beta
<b>Tested Applications</b>	ELISA,WB,IHC;WB:1:500-1:2000,IHC:1:50-1:200
<b>Form</b>	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Clonality</b>	Polyclonal
<b>Alias</b>	DKFZp686D0638; glycogen synthase kinase 3 alpha; Glycogen synthase kinase-3 alpha; GSK-3 alpha; GSK3A
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	GSK3A

## Image



Western blot analysis of GSK3 alpha/beta expression in TNF-a treated 293 whole cell lysates, The lane on the left is treated with the antigen-specific peptide.