




GSK3A Antibody

Product Code	CSB-PA205981
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49840
Immunogen	A synthesized peptide derived from Human GSK3 alpha/beta
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	GSK3 alpha/beta Antibody detects endogenous levels of total GSK3 alpha/beta
Tested Applications	ELISA,WB,IHC;WB:1:500-1:2000,IHC:1:50-1:200
Form	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	DKFZp686D0638; glycogen synthase kinase 3 alpha; Glycogen synthase kinase-3 alpha; GSK-3 alpha; GSK3A
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GSK3A

Image



Western blot analysis of GSK3 alpha/beta expression in TNF-a treated 293 whole cell lysates, The lane on the left is treated with the antigen-specific peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.