




SNAI1 Antibody

Product Code	CSB-PA162021
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95863
Immunogen	A synthesized peptide derived from Human SNAI1
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	SNAI1 Antibody detects endogenous levels of total SNAI1
Tested Applications	ELISA,WB,IHC;WB:1:500-1:2000,IHC:1:50-1:200
Form	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	SNAH, SNAI, Sna, Snail, zinc finger protein SNAI1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SNAI1

Image



Western blot analysis of SNAI1 expression in HT29 whole cell lysates, The lane on the left is treated with the antigen-specific peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.