




NFKBIB Antibody

Product Code	CSB-PA206776
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15653
Immunogen	A synthesized peptide derived from Human IkappaB-beta
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	IkappaB-beta Antibody detects endogenous levels of total IkappaB-beta
Tested Applications	ELISA,WB;WB:1:500-1:2000
Form	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	I-kappa-B-beta; IkappaBbeta; Ikb-B; Ikb-beta; IKBB; NF-kappa-B inhibitor beta; NF-kappa-BIB; NFKBIB; nuclear factor of kappa light polypeptide gene enhancer in B-cells inhibitor, beta; Thyroid receptor-interacting protein 9; TR-interacting protein 9; TRIP
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NFKBIB

Image



Western blot analysis of IkappaB-beta expression in TNF treated NIH-3T3 whole cell lysates, The lane on the left is treated with the antigen-specific peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.