



# MMP1 Antibody

<b>Product Code</b>	CSB-PA159399
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P03956
<b>Immunogen</b>	A synthesized peptide derived from Human MMP1
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Specificity</b>	MMP1 antibody detects endogenous levels of total MMP1
<b>Tested Applications</b>	ELISA, WB; WB: 1:500-1:2000
<b>Form</b>	0.01M PBS, pH 7.4 with 10 mg/ml BSA and 0.1% Sodium azide
<b>Purification Method</b>	Affinity-chromatography
<b>Clonality</b>	Polyclonal
<b>Alias</b>	22 kDa interstitial collagenase; 27 kDa interstitial collagenase; CLG; CLGN; Fibroblast collagenase; Interstitial collagenase; matrix metalloproteinase 1 (interstitial collagenase); matrix metalloprotease 1; Matrix metalloproteinase-1; MMP1; MMP1
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	MMP1

## Image

