




MMP1 Antibody

Product Code	CSB-PA159399
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03956
Immunogen	A synthesized peptide derived from Human MMP1
Raised In	Rabbit
Species Reactivity	Human
Specificity	MMP1 antibody detects endogenous levels of total MMP1
Tested Applications	ELISA, WB; WB: 1:500-1:2000
Form	0.01M PBS, pH 7.4 with 10 mg/ml BSA and 0.1% Sodium azide
Purification Method	Affinity-chromatography
Clonality	Polyclonal
Alias	22 kDa interstitial collagenase; 27 kDa interstitial collagenase; CLG; CLGN; Fibroblast collagenase; Interstitial collagenase; matrix metalloproteinase 1 (interstitial collagenase); matrix metalloprotease 1; Matrix metalloproteinase-1; MMP1; MMP1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MMP1

Image



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.