




# TFF1 Antibody

<b>Product Code</b>	CSB-PA029369
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P04155
<b>Immunogen</b>	Recombinant protein of Human TFF1
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB, IHC; WB: 1:500-1:2000, IHC: 1:50-1:200
<b>Form</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Clonality</b>	Polyclonal
<b>Alias</b>	pS2; BCEI; HPS2; HP
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	TFF1

## Image



Western blot analysis of extracts of MCF-7 cell line, using TFF1 a

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.