




TGFB1 Antibody

Product Code	CSB-PA786987
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01137
Immunogen	The antiserum was produced against synthesized peptide derived from Human TGF beta.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	TGF beta1 Antibody detects endogenous levels of total TGF beta protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:2000,IHC:1:50-1:200
Form	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol
Purification Method	Affinity-chromatography
Clonality	Polyclonal
Alias	CED; diaphyseal dysplasia 1, progressive; DPD1; TGF-beta 1; TGFB.TGF-BETA
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TGFB1

Image



Western blot analysis of extracts from HepG2 cells, using TGF beta Antibody. The lane on the left is treated with the synthesized peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.