




# GST Antibody

<b>Product Code</b>	CSB-PA187265
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	GeneID:8113787
<b>Immunogen</b>	Purified recombinant Escherichia coli BL21(DE3) GST protein fragments expressed in E.coli.
<b>Raised In</b>	Rabbit
<b>Specificity</b>	Transfected Only
<b>Tested Applications</b>	ELISA,WB;WB:1:500-1:1000
<b>Form</b>	Purified Rabbit polyclonal in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.
<b>Purification Method</b>	Affinity-chromatography
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Glutathione S Transferase; Glutathione S transferase Mu 1; Glutathione S-transferase class-mu 26 kDa isozyme; GST 26; GST; GST class mu 1; GST1; GSTM1 1; GSTM1a 1a; GSTM1b 1b; GTH4; GTM1; H B; HB subunit 4; MGC26563; MU 1; MU; SJ26 antigen; SjGST.
<b>Product Type</b>	Polyclonal Antibody

## Image



Western blot detection of GST antibody in CHO-K1 transfected with GST cell lysate and GST recombinant protein fragments?rified recombinant human GST-Lin28B & GST-Lin28A protein fragments expressed in E.coli,using GST antibody (1:10000 diluted).

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.