




ANXA5 Antibody

Product Code	CSB-PA167304
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Uniprot:P08758 GeneID:308
Immunogen	Purified recombinant Human ANXA5 protein fragments expressed in E.coli.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Monkey
Specificity	This antibody detects endogenous levels of ANXA5 and does not cross-react with related proteins.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Purified Rabbit polyclonal in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.
Purification Method	Affinity-chromatography
Clonality	Polyclonal
Alias	Anchorin CII; Annexin 5; Annexin A5; Annexin V; Annexin-5; Annexin5; AnnexinA5; AnnexinV; ANX 5; ANX A5; ANX5; ANXA5; ANXA5_HUMAN; Calphobindin I; CBP I; CBP-I; Endonexin II; ENX 2; ENX2; Lipocortin V; PAP I; PAP-I; Placental anticoagulant protein 4;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ANXA5

Image



Western blot detection of ANXA5 antibody in 3T3,HeLa,C6,COS7 and MCF7 cell lysates using ANXA5 antibody (1:2000 diluted).Predicted band size:30KDa.Observed band size:30KDa.Exposure time:20s.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.