




PHB Antibody

Product Code	CSB-PA591241
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35232
Immunogen	Peptide sequence around aa.256~260 (N-I-T-Y-L) derived from Human PHB.
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total PHB protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	Prohibitin; PHB1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PHB

Image



Western blot analysis of extracts from 3T3,HeLa,MCF,PC12 and 293 cells using PHB Antibody .