




Phospho-ABL1 (Tyr204) Antibody

Product Code	CSB-PA984556
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00519
Immunogen	Peptide sequence around phosphorylation site of tyrosine 204 (K-L-Y(p)-V-S) derived from Human c-Abl.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of Abl only when phosphorylated at tyrosine 204.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using
Clonality	Polyclonal
Alias	p150; c-ABL; EC 2.7.10.2;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ABL1

Image



Western blot analysis of extracts from COS7 cells treated with Adriamycin using Abl (Phospho-Tyr204) antibody. The lane on the right is treated with the antigen-specific peptide.

Product Modify Phospho-Tyr204