




Phospho-SRC (Ser75) Antibody

Product Code	CSB-PA928286
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12931
Immunogen	Peptide sequence around phosphorylation site of serine 75 (V-T-S(p)-P-Q) derived from Human Src.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of Src only when phosphorylated at serine 75.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using
Clonality	Polyclonal
Alias	p60-Src; c-Src; pp60c-src;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SRC

Image



Western blot analysis of extracts from COLO205 tissue using Src (Phospho-Ser75) antibody. The lane on the right is treated with the antigen-specific peptide.

Product Modify Phospho-Ser75