




Phospho-NOS1 (Ser852) Antibody

Product Code	CSB-PA197095
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29475
Immunogen	Peptide sequence around phosphorylation site of serine 852(F-N-S(p)-V-S) derived from Human n-NOS .
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of n-NOS only when phosphorylated at serine 852.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography usi
Clonality	Polyclonal
Alias	bNOS; Constitutive NOS; IHPS1; N-NOS; NC-NOS
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NOS1

Image



Western blot analysis of extracts from A549 tissue using n-NOS (Phospho-Ser852) antibody. The lane on the right is treated with the antigen-specific peptide.

Product Modify	Phospho-Ser852
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.