




Phospho-FOXO1/FOXO3/FOXO4 (Thr24/32) Antibody

Product Code	CSB-PA196080
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12778
Immunogen	Peptide sequence around phosphorylation site of threonine 24/32(S-C-T(p)-W-P) derived from Human FOXO1/3/4-pan .
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of FOXO1/3/4-pan only when phosphorylated at threonine 24/32.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography
Clonality	Polyclonal
Alias	FKHR; FOXO1; FKHR2; AFX; MLLT7
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FOXO1/FOXO3/FOXO4

Image



Western blot analysis of extracts from 293 cells treated with Serum using FOXO1/3/4-pan (Phospho-Thr24/32) Antibody. The lane on the right is treated with the antigen-specific peptide.

Product Modify	Phospho-Thr24/32
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.