




Phospho-PRKAR2B (Ser114) Antibody

Product Code	CSB-PA039825
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31323
Immunogen	Peptide sequence around phosphorylation site of Serine114(R-R-A-S(p)-V) derived from Human PKA-R2β.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Liquid
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography use
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PRKAR2B

Image



Western blot analysis of extracts from COS-7 cells treated with PMA using PKA-R2β (Phospho-Ser114) Antibody. The lane on the right is treated with the antigen-specific peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.