



Phospho-APLP2 (Tyr755) Antibody

Product Code	CSB-PA254734
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q06481
Immunogen	Peptide sequence around phosphorylation site of tyrosine 755 (P-T-Y(p)-K-Y) derived from Human APLP2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of APLP2 only when phosphorylated at tyrosine 755
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography
Clonality	Polyclonal
Alias	Amyloid protein homolog ; Amyloid-like protein 2 precursor; APP2; APPH; CDEI-box binding protein
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	APLP2
Image	 <p>Immunohistochemical analysis of paraffin-embedded human brain tissue using APLP2 (Phospho-Tyr755) antibody (left) or the same antibody preincubated with blocking peptide (right).</p>
Product Modify	Phospho-Tyr755
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.