




Phospho-CEBPA (Thr230) Antibody

Product Code	CSB-PA599746
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49715
Immunogen	Peptide sequence around phosphorylation site of threonine 230 (P-P-T(p)-P-V) derived from Human C/EBP-α.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The Antibody detects endogenous levels of C/EBP-α only when phosphorylated at Thr230.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using
Clonality	Polyclonal
Alias	CCAAT/enhancer-binding protein alpha C/EBP alpha CEBPA;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CEBPA

Image



Western blot analysis of extracts from cos-7 cells (Lane 2), using C/EBP-α (Phospho-Thr230) Antibody. The lane on the left is treated with antigen-specific peptide.

Product Modify	Phospho-Thr230
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.