




Phospho-PAK7/PAK6 (Ser602/Ser560) Antibody

Product Code	CSB-PA285604	
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.	
Uniprot No.	Q9P286	
Immunogen	Peptide sequence around phosphorylation site of Serine 602/560(R-K-S(p)-L-V) derived from Human PAK5/6.	
Raised In	Rabbit	
Species Reactivity	Human,Mouse	
Specificity	The antibody detects endogenous levels of PAK5/6 only when phosphorylated at serine 602/560.	
Tested Applications	ELISA,WB;WB:1:500-1:1000	
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.	
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography usi	
Clonality	Polyclonal	
Alias	PAK5; PAK 6; kinase PAK6;	
Product Type	Polyclonal Antibody	
Immunogen Species	Homo sapiens (Human)	
Target Names	PAK5/PAK6	
Image		Western blot analysis of extracts from HeLa cells (Lane 2), using PAK5/6 (Phospho-Ser602/Ser560) Antibody. The lane on the left is treated with antigen-specific peptide.
Product Modify	Phospho-Ser602/Ser560	
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.	