




Phospho-EPHA4 (Tyr596) Antibody

Product Code	CSB-PA959570
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54764
Immunogen	Peptide sequence around phosphorylation site of tyrosine 596 (R-T-Y(p)-V-D) derived from Human EPHA4 .
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of EPHA4 only when phosphorylated at tyrosine 596
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography usi
Clonality	Polyclonal
Alias	EK8; hEK8; EPH-like kinase 8;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EPHA4

Image



Western blot analysis of extracts from JK cells (Lane 2), using EPHA4 (Phospho-Tyr596) Antibody. The lane on the left is treated with antigen-specific peptide.

Product Modify	Phospho-Tyr596
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.