




Phospho-PHB (Thr258) Antibody

Product Code	CSB-PA097030
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35232
Immunogen	Peptide sequence around phosphorylation site of threonine 258 (N-I-T(p)-Y-L) derived from Human PHB.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of PHB only when phosphorylated at threonine 258.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography usi
Clonality	Polyclonal
Alias	Prohibitin; PHB1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PHB

Image



Western blot analysis of extracts from MCF, PC12, 293 and 3T3 cells using PHB (Phospho-Thr258) Antibody.

Product Modify	Phospho-Thr258
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.