




MMP13 Antibody

Product Code	CSB-PA981877
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P45452
Immunogen	Peptide sequence around aa.355~359 (F-W-A-L-N) derived Human from MMP-13.
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total MMP-13 protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	CLG3; Collagenase 3; MANDP1; matrix metallopeptidase 13; MMP13
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MMP13

Image



Western blot analysis of extracts from MCF,MDA231,MDA435 and HT29 cells using MMP-13 Antibody .

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.