




NTRK2 Antibody

Product Code	CSB-PA936032
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16620
Immunogen	Synthesized peptide derived from Internal of Human Trk B.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total Trk B protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	BDNF/NT-3 growth factors receptor precursor; EC 2.7.10.1; GP145-TrkB; GP145-TrkB/GP95-TrkB; kinase TrkB
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NTRK2

Image



Western blot analysis of extracts from K562 cells (Lane 2), using Trk B (Ab-705) antibody. The lane on the left is treated with synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using Trk B antibody.