





BCAR1 (Ab-410) Antibody

Product Code	CSB-PA094765
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56945
Immunogen	Synthesized non-phosphopeptide derived from Human p130 Cas around the phosphorylation site of tyrosine 410 (G-V-T-A-V).
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total p130 Cas protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	BCAR1; Breast cancer anti-estrogen resistance 1; CRK-ASSOCIATED SUBSTRATE; CRKAS;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BCAR1


Image



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using p130 Cas (Ab-410) antibody.



Western blot analysis of extracts from NIH/3T3 cells, using p130 Cas (Ab-410) antibody.



Western blot analysis of extracts from HuvEc cells using p130 Cas (Ab-410) antibody. The lane on the left is treated with synthesized peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.