





STAT5B (Ab-731) Antibody

Product Code	CSB-PA100746
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51692
Immunogen	Synthesized non-phosphopeptide derived from Human STAT5B around the phosphorylation site of serine 731(A-P-S-P-A).
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total STAT5B protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Signal transducer and activator of transcription 5B;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	STAT5B

Image



Western blot analysis of extracts from RAW264.7 cells, treated with EGF (200ng/ml 30min), using STAT5B (Ab-731) antibody.



Immunohistochemical analysis of paraffin-embedded human colon carcinoma tissue using STAT5B (Ab-731) antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.