




AKT1 (Ab-124) Antibody

Product Code	CSB-PA919481
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31749
Immunogen	Synthesized non-phosphopeptide derived from Human Akt around the phosphorylation site of serine 124 (S-G-S(p)-P-S).
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total Akt protein
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	AKT; AKT1 kinase; C-AKT; PKB; PKB-alpha
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	AKT1

Image



Western blot analysis of extracts from NIH/3T3 cells, treated with PDGF (50mg /ml, 20mins), using Akt (Ab-124) antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.