




CD19 (Ab-531) Antibody

Product Code	CSB-PA253616
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15391
Immunogen	Synthesized non-phosphopeptide derived from Human CD19 around the phosphorylation site of tyrosine 531 (D-S-Y(p)-E-N).
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total CD19 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	B-lymphocyte surface antigen B4; CD19; Leu-12;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CD19

Image



Western blot analysis of extracts from COS7 cells, treated with Serum (10%, 30mins), using CD19 (Ab-531) antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.