





FYN (Ab-530) Antibody

Product Code	CSB-PA914183
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06241
Immunogen	Synthesized non-phosphopeptide derived from Human Fyn around the phosphorylation site of tyrosine 530 (P-Q-Y(p)-Q-P).
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total Fyn protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 2.7.10.2; FYN; kinase Fyn; P59-FYN; SLK
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FYN

Image



Western blot analysis of extracts from HUVEC cells, MCF-7 cells and 293 cells, using Fyn (Ab-530) antibody.



Western blot analysis of extracts from HeLa cells (Lane 2), A549 cells (Lane 3) and HepG2 cells (Lane 4), using Fyn (Ab-530) antibody. The lane on the left is treated with synthesized peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.