




ITGB1 (Ab-789) Antibody

Product Code	CSB-PA136121
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05556
Immunogen	Synthesized non-phosphopeptide derived from Human Integrin β 1 around the phosphorylation site of threonine 789 (V-T-T(p)-V-V).
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Specificity	The antibody detects endogenous levels of total Integrin β 1 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg^{2+} and Ca^{2+}), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CD29; FNRB; Fibronectin receptor beta subunit; ITB1; Integrin VLA-4 beta subunit
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ITGB1

Image



Western blot analysis of extracts from Jurkat cells, using Integrin β 1 (Ab-789) antibody.