





JAK2 (Ab-570) Antibody

Product Code	CSB-PA281490
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60674
Immunogen	Synthesized non-phosphopeptide derived from Human JAK2 around the phosphorylation site of tyrosine 570 (G-D-Y(p)-G-Q).
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total JAK2 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 2.7.10.2; JAK-2; JAK2; Janus kinase 2; kinase Jak2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	JAK2

Image



Western blot analysis of extracts from 293 cells, treated with etoposide (25uM, 24hours), using JAK2 (Ab-570) antibody.



Immunofluorescence analysis of NIH/3T3 cells, using JAK2 (Ab-570) antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.