




ALK Antibody

Product Code	CSB-PA045830
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UM73
Immunogen	Synthesized peptide derived from Internal of Human ALK.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ALK protein.
Tested Applications	ELISA,IHC,IF;IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ALK tyrosine kinase receptor precursor; Anaplastic lymphoma kinase; CD246; EC 2.7.10.1; kinase ALK
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ALK

Image



Immunohistochemistry analysis of paraffin-embedded human brain tissue using ALK antibody.



Immunofluorescence analysis of HeLa cells, using ALK antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.