




RASGRF1 Antibody

Product Code	CSB-PA919137
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13972
Immunogen	Synthesized peptide derived from Internal of Human Ras-GRF1.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total Ras-GRF1 protein
Tested Applications	ELISA, IHC; IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CDC25; GNRP; GRF1; Guanine nucleotide releasing protein; P140 Ras-GRF
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RASGRF1

Image



Immunohistochemistry analysis of paraffin-embedded human brain tissue using Ras-GRF1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.