





# BAX (Ab-167) Antibody

<b>Product Code</b>	CSB-PA834415
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q07812
<b>Immunogen</b>	Synthesized non-phosphopeptide derived from Human Bax around the phosphorylation site of threonine 167 (F-G-T(p)-P-T).
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Specificity</b>	The antibody detects endogenous levels of total Bax protein
<b>Tested Applications</b>	ELISA, WB, IF; WB: 1:500-1:3000, IF: 1:100-1:500
<b>Form</b>	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Apoptosis regulator BAX; membrane isoform alpha; BAXA; Bcl2-associated X protein;
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	BAX

## Image



Western blot analysis of extracts from HepG2 cells, using Bax (Ab-167) antibody.



Immunofluorescence analysis of LOVO cells, using Bax (Ab-167) antibody.

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.