




PRKCA (Ab-657) Antibody

Product Code	CSB-PA589165
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17252
Immunogen	Synthesized non-phosphopeptide derived from Human PKC α around the phosphorylation site of tyrosine 657 (F-S-Y(p)-V-N).
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Specificity	The antibody detects endogenous levels of total PKC α protein
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 2.7.11.13; kinase PKC-alpha; KPCA; PKC III; PKC-A
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PRKCA

Image



Western blot analysis of extracts from MCF-7 cells and Jurkat cells, using PKC α (Ab-657) antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.