





DAPK3 (Ab-265) Antibody

Product Code	CSB-PA994602
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43293
Immunogen	Synthesized non-phosphopeptide derived from Human DAPK3 around the phosphorylation site of threonine 265 (R-M-T(p)-I-A).
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total DAPK3 protein.
Tested Applications	ELISA, WB, IHC; WB: 1:500-1:3000, IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	DAP kinase 3; DAP- like kinase; Death-associated protein kinase 3; Dlk; EC 2.7.11.1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DAPK3

Image



Western blot analysis of extracts from HuvEc cells, using DAPK3 (Ab-265) antibody.



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using DAPK3 (Ab-265) antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.