





JAK3 (Ab-785) Antibody

Product Code	CSB-PA088850
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52333
Immunogen	Synthesized non-phosphopeptide derived from Human JAK3 around the phosphorylation site of tyrosine 785 (S-D-Y(p)-E-L).
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total JAK3 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 2.7.10.2; JAK-3; JAK3; Janus kinase 3; L-JAK
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	JAK3

Image



Western blot analysis of extracts from Jurkat cells, using JAK3 (Ab-785) antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using JAK3 (Ab-785) antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.