





PAK1 (Ab-199) Antibody

Product Code	CSB-PA036079
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13153
Immunogen	Synthesized non-phosphopeptide derived from Human PAK1/2 around the phosphorylation site of serine 199 (T-K-S(p)-V-Y).
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total PAK1/2 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ADRB2; Alpha-PAK; CDC42/RAC effector kinase PAK-A; EC 2.7.11.1; P65-PAK
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PAK1

Image



Western blot analysis of extracts from K562 cells, using PAK1/2 (Ab-199) antibody.



Immunohistochemistry analysis of paraffin-embedded human placenta tissue using PAK1/2 (Ab-199) antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.