





PPP1R2 (Ab-120/121) Antibody

Product Code	CSB-PA438934
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41236
Immunogen	Synthesized non-phosphopeptide derived from Human PPP1R2 around the phosphorylation site of serine 120 and serine 121 (Q-E-S(p)-S(p)-G-E).
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total PPP1R2 protein
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	I-2; IPP-2; IPP2; Protein phosphatase inhibitor 2;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PPP1R2


Image



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using PPP1R2 (Ab-120/121) antibody.



Western blot analysis of extracts from Jurkat cells, using PPP1R2 (Ab-120/121) antibody.



Western blot analysis of extracts from JK cells (Lane 2), using PPP1R2 (Ab-120/121) antibody. The lane on the left is treated with synthesized peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.