




TGFBR1 (Ab-165) Antibody

Product Code	CSB-PA061850
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36897
Immunogen	Synthesized non-phosphopeptide derived from Human TGF β Receptor I around the phosphorylation site of serine 165 (D-P-S(p)-L-D).
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total TGF β Receptor I protein
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ADRB2; ALK-5; Activin receptor-like kinase 5; EC 2.7.11.30; ESK2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TGFBR1

Image



Western blot analysis of extracts from COLO cells and HepG2 cells, using TGF β Receptor I (Ab-165) antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.