




# TGFBR1 (Ab-165) Antibody

<b>Product Code</b>	CSB-PA061850
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P36897
<b>Immunogen</b>	Synthesized non-phosphopeptide derived from Human TGF $\beta$ Receptor I around the phosphorylation site of serine 165 (D-P-S(p)-L-D).
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Specificity</b>	The antibody detects endogenous levels of total TGF $\beta$ Receptor I protein
<b>Tested Applications</b>	ELISA, WB, IF; WB: 1:500-1:3000, IF: 1:100-1:500
<b>Form</b>	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Alias</b>	ADRB2; ALK-5; Activin receptor-like kinase 5; EC 2.7.11.30; ESK2
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	TGFBR1

## Image



Western blot analysis of extracts from COLO cells and HepG2 cells, using TGF  $\beta$  Receptor I (Ab-165) antibody.