




SNAI1 (Ab-246) Antibody

Product Code	CSB-PA041790
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95863
Immunogen	Synthesized non-phosphopeptide derived from Human SNAI1 around the phosphorylation site of serine 246 (T-F-S(p)-R-M).
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total SNAI1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	SNAH; SNAI; Sna; Snail; zinc finger protein SNAI1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SNAI1

Image



Western blot analysis of extracts from HT29 cells, using SNAI1 (Ab-246) antibody.