




CHEK1 (Ab-280) Antibody

Product Code	CSB-PA584688
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14757
Immunogen	Synthesized non-phosphopeptide derived from Human Chk1 around the phosphorylation site of Serine 280.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total Chk1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CHEK1; Chk1; kinase Chk1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CHEK1

Image



Western blot analysis of extracts from JurKat cells, using Chk1 (Ab-280) antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.