





# HDAC4 (Ab-632) Antibody

<b>Product Code</b>	CSB-PA295507
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P56524
<b>Immunogen</b>	Synthesized non-phosphopeptide derived from Human HDAC4 around the phosphorylation site of serine 632.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Specificity</b>	The antibody detects endogenous levels of total HDAC4 protein.
<b>Tested Applications</b>	ELISA,WB;WB:1:500-1:3000
<b>Form</b>	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Alias</b>	HD4; Histone deacetylase 4; KIAA0288;
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	HDAC4

## Image



Western blot analysis of extracts from HeLa/Jurkat cells, using HDAC4 (Ab-632) antibody.



Western blot analysis of extracts from HeLa cells (Lane 2), using HDAC4 (Ab-632) antibody. The lane on the left is treated with synthesized peptide.

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.