




MET (Ab-1234) Antibody

Product Code	CSB-PA955702
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08581
Immunogen	Synthesized non-phosphopeptide derived from Human Met around the phosphorylation site of tyrosine 1234 (K-E-Y(p)-Y-S).
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total Met protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	MET Hepatocyte growth factor receptor [Precursor];
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MET

Image



Western blot analysis of extracts from HepG2 cells, using Met (Ab-1234) antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.