




SNCA (Ab-125) Antibody

Product Code	CSB-PA578342
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37840
Immunogen	Synthesized non-phosphopeptide derived from Human Synuclein- α around the phosphorylation site of tyrosine 125 (E-A-Y(p)-E-M).
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total Synuclein- α protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	alpha-synuclein; NACP; Non-A beta component of AD amyloid; Non-A4 component of amyloid precursor; SYN
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SNCA

Image



Western blot analysis of extracts from rat brain and mouse brain, using Synuclein- α (Ab-125) antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.