





EGFR (Ab-1026) Antibody

Product Code	CSB-PA079221
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00533
Immunogen	Synthesized peptide derived from Internal of Human EGFR.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total EGFR protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 2.7.10.1; epidermal growth factor receptor (erythroblastic leukemia viral (v-erb-b) oncogene homolog; avian); Epidermal growth factor receptor precursor; ERBB1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EGFR

Image



Western blot analysis of extracts from HUVEC cells, using EGFR (Ab-1026) antibody.



Western blot analysis of extracts from A549 cells (Lane 2) and Hela cells (Lane 3), using EGFR (Ab-1026) antibody. The lane on the left is treated with synthesized peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.